

The Capnovet- 10 by Vetronic Services Ltd
Contents page

Introduction to the Capnovet-10.....	2
Overview and Main Working parts.....	3
Capnovet- 10 IMPORTANT NOTICE.....	4
Using the Capnovet-10 for the first time.....	5
Screen Display in Time-capnogram mode.....	5
Screen symbols, Time-Capnogram mode.....	5
Symbol descriptions – Time-Capnogram mode.....	5, 6, 7
Screen Display in Bargraph mode.....	7
Screen symbols, Bargraph mode.....	8
Symbol descriptions – Bargraph mode.....	8
The SETUP options	9, 10
Changing the SETUP options.....	10, 11
Setting the Alarms.....	12
Connecting the Capnovet-10 to a patient.....	13
Changing the screen mode.....	13
Changing the Sweep Speed.....	13
Responding to alarms.....	14
Recharging the internal batteries.....	14
Capnovet-10 recalibration.....	14
User Calibration.....	15
Factory Calibration.....	16
Appendix A.....	17

Introduction to the Capnovet-10

The Capnovet-10 is a new, microprocessor-controlled, hand-held Carbon Dioxide monitoring device designed to monitor the end-tidal CO₂ values of veterinary patients. Using side-stream technology the Capnovet-10 is suitable for a range of patients from a few hundred grammes upwards. The Capnovet-10 has two distinct modes for displaying the information. In the first mode the end-tidal value is shown as a numeric value and the real time CO₂ values are plotted as a time-capnogram. In the second mode the end-tidal value is again shown as a numeric value but the real time CO₂ values are shown as a bargraph display. In both modes the respiratory rate is also displayed in breaths per minute. The Capnovet-10 has some unique features, which are configurable by the user. The sweep speed across the screen of the time-capnogram can be altered to suit the respiratory rate of the patient being monitored. There are 5 sweep speeds, which can be changed simply and quickly by the user. Also selectable by the user is the gas withdrawal rate or sampling rate. For gas to be analysed in a side-stream device it is sucked down a small tube to the main unit. For very small animals this can start to be a problem when the gas being withdrawn approaches the tidal volume of the patient. In the Capnovet-10 this sampling rate can be altered to suit the patient. Three sampling rates can be selected and these are 50ml, 100ml and 200ml per minute. Further features of the Capnovet-10 are the ability to set alarms for high and low CO₂ settings, alarms for Apnoea time and the option to display CO₂ units in millimetres of mercury or as a percentage. With rechargeable batteries the Capnovet-10 is suitable for outdoor use or prep-room use as well as main theatre use.

Overview

Main working parts



How It Works

The Capnovet-10 continuously draws a small sample of gas from the airway attached to the patient and measures the carbon dioxide content of it by an optical means inside the unit. Gas that has been sampled is pumped out of the unit and may be returned to the anaesthetic circuit or vented in the normal way. The capnograph continuously measures the CO₂ content of the airway looking for a maximum value. When the patient starts to breathe in again the maximum found CO₂ value is shown on the screen. This is the peak end-tidal CO₂ value and will remain unchanged until the next breath is measured. At the same time, the screen also shows the instantaneous result of the measurements either as a time-based waveform (time-capnogram) or as a CO₂-based bargraph. In this way it is possible to monitor not only the end-tidal value but also the waveform during the expiratory/inspiratory phase. This waveform can provide important information on patient status. The concentration of inspired CO₂ is also shown on the display.

IMPORTANT

Because the Capnovet-10 continuously sucks air into the sensing chamber it is very important that contaminants such as dust or fluid are not allowed to enter the unit. For this reason external filters are fitted to the unit to trap fluid and particles. Under no circumstances operate the Capnovet-10 unit without these filters in place. If fluid or particles enter the sensor unit it will cease to function and require cleaning and calibration by the manufacturer. This is not covered by the guarantee.

CAPNOVET-10 IMPORTANT NOTE

Like all sidestream devices, the Capnovet-10 sucks expired gases into the detecting chamber within the instrument. For this reason sidestream devices are susceptible to moisture and other liquid contaminants and must NEVER be operated without an in-line filter. If a filter becomes blocked with repeated use or with saliva then it should be replaced. These filters can then be flushed gently with sterile saline and left to dry. Once dry they can be re-used.

The simplest way to avoid moisture contamination is to position the filter vertically above the airway so that moisture/fluid cannot drip into it.



Remember, if moisture enters the unit it will probably be necessary to return it to us for a complete recalibration. This is not covered by the guarantee and will cost approximately £200 + VAT.


Using the Capnovet-10 for the first time

Remove the Capnovet-10 from its packing and connect the sampling tubes to the main unit by mating the male luer-lock fittings with the fittings on the back of the unit. One pipe leads to the scavenging circuit and one pipe leads to the sensor unit. The sensor portion has 15mm (male & female) connectors either end to fit in line with a standard ET tube or into a face mask.

Press the ON/OFF button and a screen appears with some date information and the option to enter the SETUP menu by pressing the CHANGE button. If the CHANGE button is not pressed then the unit will wait 8 seconds (while initialisation and calibration occurs) before going straight into monitoring mode. The unit will then show a monitoring screen, either in time-capnogram mode or bargraph mode.

TIME-CAPNOGRAM MODE

Screen display in Time-Capnogram mode



In time-capnogram mode there are two elements to the screen. The main portion of the screen shows the actual capnogram trace running from left to right and the right hand portion of the screen shows the respiratory rate in breaths per minute and the end-tidal CO₂ and Inspired CO₂ values in either mmHg or %.

Screen symbols, time-capnogram mode

- | | |
|---|---|
| 1 | Battery symbol |
| 2 | Screen sweep speed |
| 3 | Respiratory rate in breaths per minute |
| 4 | End-tidal and Inspired CO ₂ value in mmHg or % |
| 5 | Free-running time capnogram |

Symbol descriptions – Time-Capnogram mode

1. Battery symbols / blocked pump



This symbol provides continuous feedback on the battery status. When the battery is fully charged the battery icon will appear filled in. As the battery discharges the battery icon will empty. The Capnovet-10 needs recharging when the battery icon appears empty as in the third icon above.



This symbol appears when the inlet or outlet to the pump is blocked and normal function is interrupted. Once the blockage is removed the normal battery symbol will return within 5-10 seconds. Readings from the unit are unreliable whilst this symbol is active.

2. Screen Sweep speed

This symbol shows the selected sweep speed for the display.

e.g. **Sweep: 2 x**

There are 5 sweep speeds selectable from 1 - 5. These equate to screen sweep speeds of 1,2,4,8 and 16 mm/second respectively. To change the screen sweep speed press the CHANGE button once, this will change to the next fastest sweep speed. Once speed 5 is reached the next key press will select speed 1 and the process will repeat. Select the sweep speed that is most consistent with the respiratory rate and visualisation of the time-capnogram.

3. Respiratory rate in breaths per minute

This value appears to the right and below the letters BPM. The respiratory rate is continuously updated with every breath. If after 15 seconds there has been no breath then the display will change – the word APNOEA will replace BPM and a reverse-video (light on dark) counter will begin counting seconds of apnoea. The reverse-video is to avoid confusion between a respiratory rate and an apnoea time.

e.g. **BPM**

23

Respiratory rate is 23
breaths per minute

APNOEA

23

Patient has not breathed
for 23 seconds

4. End tidal and Inspired CO₂ value

a) ET Value

To the right of the screen and below the CO₂ % is a value, which is the last measured end-tidal CO₂ value expressed as a percentage. If mmHg are selected as the units of end-tidal CO₂ then the value will appear under the text CO₂ mmHg.

b) IN Value

Below the End-Tidal value is another value in the same units. This is the inspired CO₂ level, which will normally read at or around zero - 0.0

5. Free-running time- capnogram.


The main area of the screen is bounded by two graticular scales between which is a graphical display of the time-capnogram. Each gradation on the scale approximates to 1% CO₂ (8mmHg). The display continuously sweeps from the left hand side to the right hand side at a rate determined by the screen sweep speed (see above). An increase in the measured CO₂ is shown as an increase in height of the time-capnogram.

Switching between Time-Capnogram mode and Bargraph mode

To switch from one screen to the other simply press the NEXT button from either screen.

BARGRAPH MODE

Screen display in Bargraph mode:



In bargraph mode the CO₂ representation is shown by a thin solid bargraph at the bottom of the screen. This bargraph is a real-time measure of the CO₂ value and will rise and fall with exhalation/inspiration. To the right of the screen at the bottom is the reported end-tidal CO₂ value for the last breath. At the top of the screen is another bargraph, which is a measure of the breaths per minute for the patient. The actual value of breaths per minute is shown to the right of this bargraph at the top right.

Screen symbols – Bargraph mode

- 1 Breaths per minute legend
- 2 Breaths per minute bargraph
- 3 Respiratory rate in breaths per minute
- 4 End-tidal and Inspired CO₂ values
- 5 Free-running CO₂ bargraph

Symbol descriptions – Bargraph mode

1. Breaths per minute legend – BPM

2. Breaths per minute bargraph

This hashed bargraph is a representation of the breathing rate of the patient. The full extent of the bargraph represents a breathing rate of 120 breaths per minute. Hence a respiratory rate of around 60 breaths per minute will appear as a bargraph of half the screen width.

3. Respiratory rate in breaths per minute

Appearing above the Breaths per Minute bargraph is a numerical value for the actual breaths per minute. After a period of apnoea for 15 seconds or the duration of the active apnoea alarm (whichever is the shorter) the apnoea time in seconds will replace the respiratory rate. This apnoea time is shown in reverse-video (light on dark) to avoid confusion with a respiratory rate.

4. End-tidal and Inspired CO₂ values

On the bottom right of the screen is a value for the last measured end-tidal CO₂ value and the Inspired CO₂ value. Depending on the units selected (see SETUP) this will either be in mmHg or as a percentage value. The End-Tidal value is updated at the end of every breath and the Inspired value at the beginning of every breath. These values will return to zero after 15 seconds of apnoea.

5. Free running CO₂ bargraph


This solid bargraph display shows the instantaneous CO₂ value on a scale from zero to full scale, where full scale is either 10% or 76mmHg.

SETUP Options

The SETUP screen

To enter SETUP mode press the CHANGE key at turn-on. If the unit is already ON you will need to turn the unit OFF and back ON again to enter SETUP mode.

The SETUP screen



SETUP options

There are two SETUP screens, a Quick Setup and an Advanced Setup.

Quick Setup

Backlight - Options: ON, OFF.

Turns the LCD screen backlight ON or OFF. Having the backlight ON will improve readability but will reduce the working battery life between re-charges.

Sound - Options: OFF, SHORT, LONG

Selects the audible beep option. Selecting OFF produces no beep on detection of a breath. SHORT gives a short low-pitched beep, LONG gives a longer higher-pitched beep.

Pump – Options: LOW, MEDIUM, HIGH

Selects the withdrawal rate of gas for sampling. Sampling rates may vary slightly but the target rates are as follows.


LOW – 50ml/minute

MEDIUM – 100ml/minute

HIGH 200 ml/minute

CO2 Unit – Options: % or mmHg

Selects the units to display CO₂ values in.



Advanced Setup

Hi ETCO₂ – Options: OFF, CO₂ values from 0.1% to 10%, or 1mmHg to 76mmHg. The value of the HIGH ET CO₂ alarm.

LowETCO₂ – Options: OFF, CO₂ values from 0.1% to 10%, or 1mmHg to 76mmHg. The value of the LOW ET CO₂ alarm.

ET Avg – The End-Tidal CO₂ value can be instantaneous or averaged over 4 or 8 breaths. 4 is the default setting and gives a smooth continuous readout

Apnoea – Options: OFF, 5 seconds to 60 seconds of apnoea time.
The allowed period of apnoea before the alarm sounds

Changing the SETUP options

To change any value, scroll down to the appropriate setting using the NEXT button and then select the required option using the CHANGE button. Once all the selections have been made scroll down to the OK option (using the NEXT button) and then press CHANGE to record all the changes made. If at any time you wish to abort the changes made, either scroll down to the CANCEL option or turn the unit OFF. To move between Quick Setup and Advanced Setup scroll to either Next or Prev and press the CHANGE button.

Setting the Backlight option

The unit may be operated with the internal screen backlight either ON or OFF. With the backlight ON the screen may be more readable in dim light. With the backlight ON the operating time of the unit is reduced as the backlight consumes appreciable battery power. In normal well-lit rooms turn the backlight OFF to save power. To change whether the backlight is ON or OFF go into SETUP mode by pressing the CHANGE button after the unit has just been turned ON. The first option indicated by the right-pointing arrow is the **Backlight** option. Pressing the CHANGE button now will toggle the setting between ON and OFF. Once the selection has been made, scroll down to OK at the bottom of the screen and press CHANGE to exit and save the changes.

Changing the Sound option

During patient monitoring each breath can be marked by a beep, which can be a high or low pitch or the beeps may be turned OFF. To change the respiratory beep settings go into SETUP mode by pressing the CHANGE button after the unit has just been turned ON. Use the NEXT button to scroll down to '**Sound**'. Pressing the CHANGE button will iterate through the options of OFF, SHORT or LONG. Choose the required option and then scroll down to OK at the bottom of the screen and press CHANGE to exit and save the changes.

Changing the Pump option

The Capnovert-10 can be set up to withdraw 50ml, 100ml or 200ml per minute from the patient's airway. In any side-stream system there is a trade-off between sampling rate and signal quality. With more gas presented to the sensor the better the signal quality will be. However for smaller patients the higher withdrawal rates may represent samples that start to approach the tidal volume of the patient itself. In these circumstances gas in the expiratory limb will be pulled into the unit mixing with the exhaled gas and giving a falsely depressed reading. For these patients a low sampling rate will give a better end-tidal value at the expense of a slightly rounded waveform.

To change the pump value go into SETUP mode by pressing the CHANGE button just after the unit has been turned ON. Scroll down to the **Pump** option using the NEXT button. Pressing the CHANGE button will now iterate between LOW, MEDIUM and HIGH, representing the 50, 100 and 200 ml per minute sampling rates respectively. When the rate has been selected, scroll down to OK using the NEXT button and press CHANGE to exit and save the changes.

Changing the CO₂ units

The value for the end-tidal CO₂ can be expressed either as a percentage or in millimetres of mercury (see Appendix A for conversions). To change the CO₂ units, go into SETUP mode by pressing the CHANGE button just after the unit has been turned ON. Scroll down to the **CO₂** option and press the NEXT button to iterate between % and mm Hg. When the selection has been made scroll down to OK using the NEXT button and press CHANGE to exit and save the changes.

Please note that when changing between CO₂ units the alarm values for mmHg and %CO₂ are independent of each other and need to be set explicitly for the units selected.

Setting the Alarms

TO set an alarm, turn the unit ON and immediately press the CHANGE button. The SETUP screen appears. Using the NEXT button scroll down to the alarm parameter that you want to change.

Setting an upper limit for ET CO₂

An alarm can be set to sound if the patient's CO₂ value exceeds the value in the HI ET CO₂ setting. Scroll down to **HI ET CO₂** using the NEXT button. Then use the CHANGE button to alter the value. Note that you cannot set a HI ET CO₂ value lower than the LO ET CO₂ value. Pressing the CHANGE button once increments the value by 0.1% or 1mmHg. If the CHANGE button is pressed and held, then the numbers increment automatically. When the numbers reach the upper limit the values pause at the OFF state for a couple of seconds before incrementing again.

Setting a lower limit for ET CO₂

An alarm can be set to sound if the patient's CO₂ value falls below the value in the LO ET CO₂ setting. Scroll down to **LO ET CO₂** using the NEXT button. Then use the CHANGE button to alter the value. Note that you cannot set a LO ET CO₂ value higher than the HI ET CO₂ value. Pressing the CHANGE button once increments the value by 0.1% or 1mmHg. If the CHANGE button is pressed and held, then the numbers increment automatically. When the numbers reach the upper limit the values pause at the OFF state for a couple of seconds before incrementing again.

Please note that when changing between CO₂ units the alarm values for mmHg and %CO₂ are independent of each other and need to be set explicitly for the units selected.

Setting the Averaging for End-Tidal values

The End-Tidal value can be displayed as an instantaneous value (breath to breath) or as an averaged value over 4 or 8 breaths. This option is selected by using the NEXT button to scroll down to **Next** in the Quick Setup screen and then pressing the CHANGE button. This takes you to the Advanced Setup screen. Scroll down to the **ET avg** line and press CHANGE. The value changes between 1, 4 and 8 repeatedly. Once the required selection has been made, scroll down to **Ok** on the bottom line and press CHANGE.

Setting an APNOEA alarm time

An alarm can be set to sound if the patient fails to take a breath. The apnoea alarm time limits are between 5 and 60 seconds in one-second increments. Scroll down to **Apnoea** using the NEXT button. Then use the CHANGE button to alter the value. To turn the function OFF again or to scroll to a lower value than shown, then press and hold the CHANGE button. The values will auto-increment until the upper limit is reached when they will pause on the OFF value before incrementing again.

Connecting the Capnovet-10 to a patient

The Capnovet-10 comes complete with an in-line 15mm connector that will connect directly to the 15mm connector on the endo-tracheal tube or to the back of a standard face mask. The Capnovet-10 uses side-stream technology to determine the patient's CO₂ value. This means that a continuous sample of gas is removed from the patient's airway and taken to the main unit via the small-bore tubing. Once sampled, this gas can be dealt with in two ways. It may be vented to the normal scavenging system by connecting the outlet port to the scavenging system or it may be returned to the inspiratory side of the circuit in a rebreathing system. It is important that this sampled gas is dealt with correctly; otherwise it may effectively act as a source of gas contamination to the environment. Also supplied with the Capnovet-10 is a set of ET tube connectors with side port connectors allowing the sampling tube to be connected at the very end of the ET tube. This significantly improves the waveform seen in the smaller animals. With the sampling connector in place turn the Capnovet-10 ON, make any necessary changes to the SETUP options and proceed with monitoring the patient.

Changing the screen mode

Every time the Capnovet-10 is turned OFF, the display mode you were viewing is remembered and will re-appear when you next turn it ON. To change the display mode press the NEXT button. This will allow you to toggle back and forth between display modes. Note that when in Bargraph mode the CHANGE button will have no effect.

Changing the Sweep Speed

With the Time-capnogram mode display the rate at which the display is swept across the screen can be altered. Simply pressing the CHANGE button whilst in this mode will advance the sweep rate. Five sweep rates are available from 1-5.

Responding to alarms

When an alarm condition is detected (Apnoea time exceeded or end-tidal CO₂ value detected outside of the alarm limits) an intermittent alarm will sound. If the alarm condition is self-righted, as for example the patient taking a breath during an Apnoea alarm, the alarm will stop. Otherwise the alarm will continue to sound until the condition has been righted or the alarm is temporarily silenced. The Capnovet-10 has an in-built 60-second alarm silence function, which is activated during an alarm condition. While an alarm is sounding pressing either the NEXT or CHANGE button will silence the alarm for 60 seconds only. At the end of this 60-second period the alarm will again sound if the alarm condition has not been rectified. Note also that Apnoea takes priority in the alarms. If an alarm has been silenced for a HI ET CO₂ cause and the patient is then apnoeic, the apnoea alarm will then sound even though 60 seconds may not have elapsed since the previous alarm was silenced.

Recharging the internal batteries

The Capnovet-10 uses internal NimH batteries which maybe re-charged using the supplied charger. NimH batteries offer a high energy density power source with quick re-charging times and no memory effect, as seen with NiCad batteries. When the battery symbol shows that the battery levels are low (see page 5, Battery Symbols) plug the charger into the mains and connect the small plug to the unit using the connector on the right-hand side of the unit. The time taken to fully re-charge the batteries is around 5 hours but will depend on how run down they are. The unit can still be used whilst the batteries are charging although the re-charging time will be extended.

The Capnovet-10 has an in-built intelligent charging circuit so that even if the charger is left connected for extended periods no damage will result.

Be aware that if the batteries become completely exhausted then the user settings will be lost. These settings include the alarm limits, pump speeds, sound and backlight options. Also if the unit has been re-calibrated the user calibration settings will also be lost.

Capnovet-10 recalibration

The Capnovet-10 has been factory calibrated prior to despatch. To maintain the accuracy of your Capnovet-10 it is important that the unit be re-calibrated once a year. On the opening screen is a re-calibration date, which is when the next calibration is due. Calibration can be performed in two ways: User calibration and factory calibration.


User calibration

User calibration requires the use of a 5% Carbon Dioxide mixture and access to fresh air. 5% Carbon Dioxide mixtures can be obtained from a variety of sources. If you do not have access to a calibrated 5% source, then contact the supplier of your Capnovet-10 for information on how to obtain some.

Procedure for User calibration

Hold down the CHANGE button firmly whilst turning the unit ON. A screen will appear which looks like the screen on the next page:

User Calibration main screen



Initially the cursor arrow will be against the **Factory reset option**. Press CHANGE to enter this option. The current user delta values will be shown as 3 values on the screen. If the user deltas have not been set before then they will read

1 1 0

Press CHANGE to reset the delta values to 1,1,0 and then press NEXT to exit this screen.

Now move the cursor, using the NEXT button to the **Air adjustment** option and press CHANGE. Make sure the unit is sampling from clean fresh air outside of the building. Room air can become stale and CO2 levels can rise to 0.4 or 0.5 %, which is too high for accurate re-calibration. When the screen value is steady, press the CHANGE button. The word "Calibration" will appear and the running dots after it indicate the unit is busy. After a few seconds the screen display will change to 0.0% and the zero point is set. Press NEXT to leave this option. Now move the cursor using the NEXT button to the **5% CO2 adjustment** option and press CHANGE. You now need to flow your 5% carbon dioxide gas through the ET connector/sensor at a rate of around 200mL/minute. Wait for the screen value to stabilise and then press the CHANGE button. The same process of

auto-calibration will be performed as is done with the Air calibration. Wait for the screen value to change to 5.0% and then press the NEXT button to leave this option.

Finally, leave this User Calibration facility by placing the cursor next to the **Done** option and pressing the CHANGE button.

Calibration is complete.

Factory Calibration

This requires sending the Capnovet-10 unit back to your supplier. Contact the supplier listed on the back of this manual to arrange this service.

Appendix A

Specifications & conversions

The Capnovet-10 can display end-tidal CO₂ values in either mmHg (millimetres of mercury) or as a percentage (%).

To convert from one set of units to the other use the following:

1% CO₂ is equivalent to 7.6mmHg

1mmHg is equivalent to 0.132% CO₂

(Note this is based on conditions at STP of 760mmHg per atmosphere)

Technical Specifications

Mechanical

Dimensions	25cm x 12cm x 10cm
Weight	625g
Construction	Rugged Polyamide construction in dark grey

General

Screen size	60mm x 38 mm viewing area
Screen type	Reflective LCD 128 x 64 pixels with integral backlight
Sensor connection	3mm x 60cm nylon tubing connecting to 15mm T-connector
Sampling rate	20Hz

Electrical Supply	Internal rechargeable NimH batteries
Power	1800mW
Current consumption	300mA nominal
Battery life	6 hours continuous use
Recharging time	5 hours

CO ₂ detection range	0-10% 0-76mmHg
---------------------------------	-------------------

Alarms	High/Low CO ₂ levels, apnoea time
Respiratory rate range	1-120 breaths per minute
Apnoea alarm settings	5-60 seconds

Display modes	1) Time-capnograph, capnometer & respiratory rate 2) CO ₂ bargraph, capnometer & respiratory rate
Sweep speeds	1,2,4,8,16 mm/second screen sweep speed

Side-stream withdrawal rates	50, 100 & 200 ml/minute for LOW, MEDIUM and HIGH respectively
------------------------------	---

Extras for your Capnovet- 10

Part Number: VSL-1948 Description: 15M-15F Swivel stem 6mm	Part Number: VSL - 2710 Description: Leur Elbow clip for 7.6mm Port – latex free	Part Number: VSL - 2731 Description: Monitor Line 1.8M Male-Female Leur Lock
 VSL-1948	 VSL - 2710	 VSL - 2731

Part Number: VSLETSP + Size Description: ET Connector with side port available in sizes 2.5, 3.0 or 3.5mm Male end	Part Number: VSLETSP + Size Description: ET Connector with side port available in sizes 4.0 or 5.0 mm Male end	Part Number: VLSRP15 Description: Filter for Capnovet lines
 VSLETSP2.5 VSLETSP3.0 VSLETSP3.5	 VSLETSP4.0 VSLETSP5.0	 VLSRP15

For other products and consumables see our web page:
<http://www.vetronic.co.uk>