



Reference Notes for Capnography

Normal Values

Normal End-tidal CO₂ values: 35-40mmHg (4.5%- 5.5%)

Normal Inspired CO₂ value 0.0% (0mmHg)

Suggested Alarm Limit Settings:

High End-Tidal CO₂ alarm setting 6.5% (50mmHg)

Low End-Tidal CO₂ alarm setting 2.0% (15mmHg)

Inspired CO₂ alarm setting 1.0% (8mmHg)

Corrective Action:

High end-Tidal CO₂ levels:

Check condition and presence of Soda Lime. Replace if necessary

Check for sufficient fresh gas flow from anaesthetic machine

Give a few extra manual breaths to increase the Minute Volume of the patient and reduce the end-tidal CO₂ value

Low End-Tidal CO₂ levels:

Most commonly occur during mechanical ventilation due to hyper ventilation. In these cases reduce the ventilation Minute Volume – either a slower respiratory rate and/or a reduced tidal volume

Can occur due to tachypnoea following induction or light anaesthesia.

Unless prolonged, no cause for alarm. Check patient depth of anaesthesia.

High Inspired CO₂ value

Exhausted or absent soda lime. Check and replace if necessary.

Insufficient fresh gas flow, leading to re-breathing. Increase FGF

Blocked/sticky expiratory valves (especially on ADE systems) leading to rebreathing.

Low Flow Anaesthesia

During Low Flow Anaesthesia with a closed system the oxygen delivery must at least meet the metabolic requirement for the animal.

Oxygen requirement:

Basic metabolism = 10mls/kg/minute.

(Varies with body temperature and metabolic rate)

e.g. 30kg Labrador needs a minimum of 300mls/minute oxygen flow.

In practice you would set this slightly higher than this at say 4-500mls/min. During low flow anaesthesia monitor a) breathing circuit oxygen (and patient SpO₂) levels for potentially hypoxic mixtures and b) circuit CO₂ levels indicating absorbent exhaustion.