

*Quality Anaesthetic monitoring
and Ventilation equipment*



Capnography - Principles & Pitfalls

**Keith Simpson BVSc MRCVS
MIET(Electronics)
Abbotskerswell, Devon.**

*Designed
by Vets, for Vets*



What is Capnography?

- Capnography is the science of monitoring carbon-dioxide levels in the respiratory gases of an animal
- Capno = Greek for smoke/vapour
- The equipment used is a Capnograph

*Designed
by Vets, for Vets*



Capnography – Sidestream or Mainstream?

- Mainstream Unit – a device that samples the CO₂ levels in-line. There is no delay in the capnogram trace. Gives a fast response. Fixed volume of dead space
- Sidestream Unit – a device that extracts a sample of the airway gas and performs the analysis inside the machine. Can be very small dead space. Needs a low sampling rate for tiny animals.

*Designed
by Vets, for Vets*




Capnography - Sidestream



*Designed
by Vets, for Vets*



Capnography - Mainstream



IRMA mainstream CO2 Monitor

*Designed
by Vets, for Vets*



Sidestream/ Mainstream

- Demonstration

*Designed
by Vets, for Vets*



Why look at CO₂?

- Measurement is a relatively non-invasive technique and depending on the method used can be performed with conscious animals
- Exhaled carbon dioxide levels give an easy way of assessing the state of the respiratory and or cardiovascular system
- Capnography is probably the single most useful monitoring modality
- Virtually no false errors – high reliability

*Designed
by Vets, for Vets*



What information do we get?

- Information from Capnography can be broken down into 4 levels, each with increasing degrees of information

*Designed
by Vets, for Vets*



Capnography - Information

- **Level 1:**

- Breathing or not, i.e. apnoea monitor
- Respiratory rate

*Designed
by Vets, for Vets*



Capnography - Information

- **Level 2:**

- Expired CO₂ levels
- Inspired CO₂ levels

- From these parameters we can now begin to deduce the state of the patient with regard to respiration i.e. normocapnic, hypocapnic or hypercapnic

*Designed
by Vets, for Vets*



Capnography - Information

- **Level 3:**

- Waveform profile
- There are 4 recognised parts to a typical capnogram, each one having characteristics that impart specific information

*Designed
by Vets, for Vets*



Capnography - Information

- So there is a huge amount of information available from a standard capnograph.
- In viewing the wealth of information from a capnograph we should not lose sight of the most basic information..

*Designed
by Vets, for Vets*



Capnography - Principles

- Boils down to the fundamental fact that respiring cells consume oxygen and produce carbon dioxide as a by-product and that this carbon dioxide is removed by the lungs
- This continuous production of carbon dioxide is the essence of determining patient status
- The summary of the above is that all situations can be reduced to a balance of production and removal of CO₂

*Designed
by Vets, for Vets*



Capnography - Principles

- Essentially for a resting animal the rate of production of CO₂ is constant. It is varied by the metabolic rate of the animal, so in some circumstances can change.
- Hyperthermia will lead to an increase in CO₂ production.
- Conversely hypothermia will lead to a reduction in CO₂ production

*Designed
by Vets, for Vets*




Capnography - Visualisation

- Consider the production of CO₂ as an open water tap filling a large vessel
- In this concept the vessel is the lungs and the open tap is the constant production of CO₂
- Consider the removal of CO₂ as a man removing water from the vessel with a bucket

*Designed
by Vets, for Vets*



Capnography - Visualisation



Rate of removal = rate of delivery NORMOCAPNOEA

*Designed
by Vets, for Vets*



Capnography - Visualisation


•In the previous slide:

- The running tap represents a constant delivery of CO₂
- The big bucket represents the level of CO₂ in the lungs
- The men represent the rate of removal of CO₂ from the lungs

*Designed
by Vets, for Vets*



Capnography - Visualisation



Rate of removal = rate of delivery NORMOCAPNOEA

*Designed
by Vets, for Vets*



Capnography – the science

- Therefore to keep an animal normocapnic we must ensure that the rate of breathing is sufficient to remove the CO₂ delivered to the lungs
- We will return to this concept later

*Designed
by Vets, for Vets*



Capnography

- Now that we understand where the information is coming from, what can we do with this information, and what does it mean for the animal?

*Designed
by Vets, for Vets*



Capnography - Level 2 information

- Normal values
- Inspired CO₂ : 0%
- Expired CO₂: 4.5% (35 mmHg)
- The following slides assume a zero level for inspired CO₂. We will come to the effects of increased inspired CO₂ later

*Designed
by Vets, for Vets*




Capnography - Hypercapnia

- Increased levels of end-tidal (alveolar) CO₂.
 - 5-6% begin to question why
 - > 6% take some action
- This means that the level of CO₂ in the lungs after equilibration with pulmonary blood is too high, which itself means that the removal of CO₂ is insufficient
- This usually means that the animal is under-ventilated. This may indicate respiratory depression or some other cause (lung pathology)
- Rarely does it reflect an increase in metabolic CO₂ production

*Designed
by Vets, for Vets*

Capnography - Hypercapnia



Rate of removal < rate of delivery: Hypercapnia

*Designed
by Vets, for Vets*



Capnography - Hypercapnia

- What's wrong with hypercapnia?
- Indicates hypo-ventilation, so reduced oxygen exchange, and anaesthetic exchange
- Cerebral vasodilation increasing cerebral blood flow and intracranial pressure
- Peripheral vasodilation by direct effect on vessels
- Cellular acidosis develops rapidly affecting cellular function
- Central depression at very high levels of $p\text{CO}_2$
- Stimulation of the sympathetic nervous system resulting in tachycardia, peripheral vasoconstriction and sweating

*Designed
by Vets, for Vets*



Capnography - Hypercapnia

- In a clinical situation this indicates a patient that has insufficient minute volume ventilation.
- It can indicate a patient is too deep
- This patient needs an increase in minute volume ventilation: Give supplementary breaths by hand by squeezing the bag or begin mechanical ventilation

*Designed
by Vets, for Vets*



Capnography - Hypercapnia

- Decreased levels of end-tidal (alveolar) CO₂.
- This means that the level of CO₂ in the lungs after equilibration with pulmonary blood is too low, which itself means either that the removal of CO₂ is excessive or that CO₂ delivery is reduced

*Designed
by Vets, for Vets*




Capnography - Hypercapnia

•Possible causes:

- Excessive ventilation
- Reduced CO₂ production – hypothermia
- Reduced delivery of CO₂ to the lungs

*Designed
by Vets, for Vets*

Capnography - Hypercapnia




Rate of removal > rate of delivery: Hypocapnia

*Designed
by Vets, for Vets*



Capnography - Hypercapnia



Reduced delivery of CO₂: Hypocapnia

*Designed
by Vets, for Vets*



Capnography - Hypercapnia

- **Reduced delivery of CO₂**
 - Can be either reduced cardiac output or decreased production
 - Differentiate between the two by degree of change. With a change in cardiac output the effect on the capnogram is dramatic

*Designed
by Vets, for Vets*



Capnography - Hypercapnia

- **What's wrong with Hypocapnia?**
- Inhibition of respiratory drive
- Shift of the haemoglobin oxygen dissociation curve to the left (impairing peripheral oxygen unloading)
- Depression of myocardial contractility
- Cardiac arrhythmias
- Less of an evil than Hypercapnia

*Designed
by Vets, for Vets*



Capnography – Waveforms

- So far all the information discussed has been available from a capnometer, i.e. the information has been purely numerical.
- We will now look at the information imparted by the capnogram itself

*Designed
by Vets, for Vets*




What can a Capnogram tell us?

- In addition to the End-Tidal value the Capnogram shows the levels of CO₂ during the complete respiratory cycle
- The Capnogram waveforms provide a lot of information

*Designed
by Vets, for Vets*




Typical Capnogram



*Designed
by Vets, for Vets*



What do the different phases represent?




Phase 0 – Inspiration. Rapid fall in CO₂ levels as new breath taken

*Designed
by Vets, for Vets*



Phase I




Phase I – Beginning of expiration and the elimination of gas from anatomical dead space

*Designed
by Vets, for Vets*



Phase II




Phase II – Alveolar gas mixes with dead space

*Designed
by Vets, for Vets*



Phase III




Phase III – Main component of expiration. Has a typical positive slope

*Designed
by Vets, for Vets*



Normal Capnogram



Expired CO₂: 5% 38mmHg
Inspired CO₂: 0% 0 mmHg
Exp – Insp: 5% 38 mmHg

*Designed
by Vets, for Vets*




Increasing End-tidal CO₂

- CO₂ delivery to the lungs exceeds removal
- Need to increase the minute volume to “blow off” more CO₂
- If animal is not on a ventilator, then start IPPV
- If animal is on a ventilator then increase the Respiratory Rate
- Could increase the Tidal Volume but that might change the degree of lung inflation
- Why has this occurred? Was Minute Volume too low or is there another factor? Hyperthermia, changes in Cardiac Output

*Designed
by Vets, for Vets*



Capnogram 2



Expired CO₂: 7.5% 55 mmHg
Inspired CO₂: 0% 0 mmHg
Exp – Insp: 7.5% 55 mmHg

*Designed
by Vets, for Vets*




Falling End-Tidal CO₂

- CO₂ removal rate exceeds the delivery rate to the lungs
- Need to reduce the minute volume
- If the animal is not on a ventilator then it is probably hyperventilating. Possible response to pain - may need increased depth of anaesthesia
- If the animal is on a ventilator then the minute volume is too high. Reduce the minute volume by decreasing the respiratory rate.
- Another cause could be....
- Fall in cardiac output so delivery of CO₂ to the lungs is decreased. Causes – compression of pulmonary artery by surgeon, PEEP, shock
- Hypothermia. Results in reduced CO₂ production

***Designed
by Vets, for Vets***



Capnogram 3



Expired CO₂: 3% 25mmHg
Inspired CO₂: 0% 0 mmHg
Exp – Insp: 3% 25 mmHg

*Designed
by Vets, for Vets*



Dramatic Drop in CO₂ level between 2 breaths

- Such a sudden fall in CO₂ level must be associated with a fall in cardiac output. This can happen when the surgeon “embarrasses” heart function or occludes the pulmonary artery.
- Something to look out for in ruptured diaphragm repair

*Designed
by Vets, for Vets*



Capnogram 4




Expired CO₂: 6% 45 mmHg then 3% 25 mmHg
Inspired CO₂: 0% 0 mmHg

*Designed
by Vets, for Vets*



Capnogram 5



Expired CO₂: 7% 55 mmHg
Inspired CO₂: 2% 17 mmHg
Exp - Insp: 5% 38 mmHg

*Designed
by Vets, for Vets*




Rebreathing

- Live Demo on a standard Y-piece..

*Designed
by Vets, for Vets*



The T-piece and gas flow



The FGF must be sufficient to meet the needs of the patient during the inspiratory phase, otherwise rebreathing will occur.

It is very easy to under-estimate the FGF requirement.

How do you calculate the required FGF?

*Designed
by Vets, for Vets*



The T-piece and gas flow

- To calculate the FGF needed for inspiration:
- Calculate the volume required in 1 second and then x 60
- e.g. 5kg cat. T.V. 50mls. Inspiratory time (IT) 1 second = $50\text{mls/second} \times 60 = 3000\text{mls/minute} = 3.0\text{L/minute}$
- Note that it IT is actually 0.8s: 3.75L/min

*Designed
by Vets, for Vets*



The T-piece and gas flow

- If you have a Capnograph you can use it to set the required FGF.
- Note that the FGF will vary over the course of a procedure because tidal volume will change

*Designed
by Vets, for Vets*



Beware hand ventilation

- You must be aware of the type of anaesthetic circuit in use when hand ventilating.
- Two basic types:

*Designed
by Vets, for Vets*



Hand Ventilation

- Re-breathing bag in supply limb e.g. Magill or Lack. High risk of re-breathing when performing IPPV
- Re-breathing bag in waste limb e.g. T-piece or Bain (modified Bain). Less risk of re-breathing during IPPV.

*Designed
by Vets, for Vets*




Hand Ventilation

- Check the Capnograph profile during IPPV and increase FGF if there is any sign of re-breathing.

*Designed
by Vets, for Vets*




Delta Capnography – Nonin 9874V



***Designed
by Vets, for Vets***



Capnogram 6



*Designed
by Vets, for Vets*



Skewed Capnogram

- The Capnogram is skewed.
- The total cycle time is unchanged but the time for phase II is extended
- Phase II is the expiration of gas from the lungs, mixing with dead space gas.
- Elongation of phase II is most commonly caused by obstruction in the expiratory circuit

*Designed
by Vets, for Vets*



Spiked Capnogram



*Designed
by Vets, for Vets*



Spiked Capnogram

- The spike represents an isolated pocket of high percentage CO₂ gas during the expiratory period.
- Known as gas trapping this can be due to:
 - Single-bronchus intubation
 - Mucus plugs

*Designed
by Vets, for Vets*



Capnography Summary

- Easy to use
- Lots of information available from basic through to detailed
- The single most informative piece of equipment for anaesthetised animals
- Remember to consider production versus removal in all cases where CO₂ output changes

*Designed
by Vets, for Vets*



Thank-you for Listening

If you have any further issues or enquiries about Vetronic Services LTD and their products please do not hesitate to contact us on the below contact details:

Tel: 01626 365505

Fax: 0870 129 4705

Email: enquiries@vetronic.co.uk

Web: www.vetronic.co.uk



*Designed
by Vets, for Vets*