

*Quality Anaesthetic monitoring
and Ventilation equipment*



ECG Monitoring

**Keith Simpson BVSc MRCVS
MIET(Electronics)
Abbotskerswell, Devon.**

*Designed
by Vets, for Vets*



Getting a good clean ECG trace

- The Equipment
- Methods of attachment
- Assessment

*Designed
by Vets, for Vets*



ECG cable selection

- Shielded or unshielded
- Makes a lot of difference to the ability to obtain a good ECG
- Many medical ECG cables are shielded only as far as the split. After that the cable is unshielded
- Susceptibility to mains interference is greatly increased.

*Designed
by Vets, for Vets*



UK lead convention

- RED – Right Fore
- YELLOW – Left Fore
- GREEN – Left Hind
- BLACK – Right Hind (reference electrode)

*Designed
by Vets, for Vets*



USA lead convention

- RED – Left Hind
- WHITE – Right Fore
- BLACK – Left Fore
- GREEN – Right Hind (reference electrode)

*Designed
by Vets, for Vets*



Methods of attachment

- **Primary aims**

- Maximise the surface area of attachment
- Minimize the contact resistance
- Minimize noise pickup

*Designed
by Vets, for Vets*



Maximize the surface area of attachment

•Options


- Crocodile clips
- Sticky ECG electrodes
- Limb Plate Electrode
- Hypodermic needles

*Designed
by Vets, for Vets*



Crocodile Clips

- **Can be painful**
 - Use modified clips
 - 'Bite' as much skin as possible
 - Spreads the load
 - Reduces pain
 - Maximises the contact
- Need a conducting medium
- Much better with shaved area



*Designed
by Vets, for Vets*



Sticky ECG electrodes

- Variety of sizes but the contact area is essentially the same
- Need well prepared dry skin
 - Clip area larger than ECG pad
 - De-grease with surgical spirit
 - Allow to dry and then apply pad
- Takes time to apply so good for long ops or ambulatory monitoring
- Can get reasonable results from the pads of the feet but skin resistance is high so subject to noise and movement artefact
- Good preparation is imperative otherwise get mediocre results

*Designed
by Vets, for Vets*



ECG pad types



*Designed
by Vets, for Vets*



Limb plate electrodes

- Very easy to apply
- Non painful
- Large surface area so good signal pickup
- No need to clip hair
- Good long term results (hours)



*Designed
by Vets, for Vets*



Hypodermic Needles

- Fairly well tolerated – 21g or 23g
- Bypasses skin resistance but small contact area so can produce unreliable results, particularly with movement
- Good for e.g. birds where skin is v fragile and cannot use pads or LPE
- Push needle under skin and out again, then attach crocodile clip to free end – picture.....

*Designed
by Vets, for Vets*



Minimize contact resistance

- Hair provides distinct barrier to contact, even when soaked in spirit
- Remove hair – pluck or shave
 - Anaesthetised dogs/cats: Pluck hair between pads and use modified crocodile clip for quick connection
- For crocodile clips and LPE apply contact medium *before* attaching electrode
- Use a conducting medium
 - Surgical spirit/isopropyl alcohol
 - Electrolyte solution/Saline
 - Cardiac gel/Ultrasound gel (remove after use)

***Designed
by Vets, for Vets***

<i>Agent</i>	<i>Pros</i>	<i>Cons</i>	<i>Contact Time</i>
Surgical Spirit	Promotes good contact. Cheap. Easily available	Smell. Irritant. Short contact time	10-15 minutes
Isopropyl Alcohol	Promotes good contact. Reasonably cheap	Smell. Irritant	30-40 minutes
Electrolyte Solution, e.g. Signa Spray by Parker	Non-irritant. Good contact. Low/zero odour	Less readily available. Increased cost	1 hour+
Cardiac Gel	Good sustained contact	Less easily available. Increased cost. Messy. Clogs up electrodes.	2 hours+
Self adhesive Electrodes	Clean. New contact every time. Simple to use	More expensive. Must prepare area. Not re-usable	2 hours+ Up to days in some cases



Minimize contact resistance

Don't even try with a set of clips like this



*Designed
by Vets, for Vets*



Minimize Noise Pickup

- So now you have a well attached ECG set of cables. Is that all you need for a good ECG?

- **NO!**

- Must reduce noise pickup as well

*Designed
by Vets, for Vets*




Cable Position

- Cable position is extremely important
- Due to phenomenon of Common Mode Rejection
- Simply put – For best rejection of noise all leads must “see” the same interference and so follow the same path

*Designed
by Vets, for Vets*



Common Mode Rejection



*Designed
by Vets, for Vets*



Improving Common Mode Rejection

- Minimise contact resistance
- All leads should be same length especially if unshielded
- All leads should run as close together as possible




*Designed
by Vets, for Vets*



Effect of cable spread

- The two pictures show the ECG of the patient with spread cables and with parallel cables



Demonstration of Effect of Cable Spread




*Designed
by Vets, for Vets*



So you can not get
anything but noise

What do you do?

- Test the machine
- Remove Causes of interference



*Designed
by Vets, for Vets*



Test the Machine


- Two very simple tests to assess the ECG machine
 - All Leads together
 - Individual Lead test

*Designed
by Vets, for Vets*



All leads together

- On all Lead settings of the machine there should be a completely flat line. If not, there is a fault




*Designed
by Vets, for Vets*



Individual Lead Tests

- Hold all four leads in one hand and fan out the clips



*Designed
by Vets, for Vets*



Individual Lead Tests

- Touching any one of these leads (except Right Hind) should cause massive stylus movement or screen trace deflection.
- Touch each in turn. A properly functioning machine will have the following response

*Designed
by Vets, for Vets*



Individual Lead Tests

•Lead I

- Response with Right Fore and Left Fore
- No response with Left Hind or Right Hind

•Lead II

- Response with Right Fore and Left Hind
- No response with Left Fore or Right Hind

•Lead III


- Response with Left Fore and Left Hind
- No response with Right Fore or Right hind

*Designed
by Vets, for Vets*




Deducing a lead problem from the ECG result


Right Hind contact problem



Lead I



Lead II




Lead III




Deducing a lead problem from the ECG result


Right Fore contact problem



Lead I



Lead II




Lead III (Does not involve RF)



Deducing a lead problem from the ECG result


Left Fore contact problem



Lead I



Lead II (Does not involve LF)




Lead III




Deducing a lead problem from the ECG result


Left Hind contact problem



Lead I (Does not involve LH)



Lead II



Lead III



Sources of Interference

- Mains induced noise
- Movement
- Mobile phones
- Newer digital Cordless phones

*Designed
by Vets, for Vets*



Thank-you for Listening

If you have any further issues or enquiries about Vetronic Services LTD and their products please do not hesitate to contact us on the below contact details:

Tel: 01626 365505

Fax: 0870129 4705

Email: enquiries@vetronic.co.uk

Web: www.vetronic.co.uk



***Designed
by Vets, for Vets***