

Changes in the supply of Nomolines for Vetronic Services Equipment

The supply of the traditional Nomoline (108210) terminating in a male luer connector will be discontinued at the end of June 2020

This is a manufacturing change imposed by the supplier and manufacturer to improve patient safety by preventing incorrect connection of supply and sampling lines to patients.

For our veterinary customers there will be 3 possible alternatives:

Order code 3827

This is the Adult sampling line which terminates in a 15mm Male-15mm Female straight connector and has a full 15mm central bore. This is suitable for patients of 5kg and above.

Price: £12.50 each



Order code 3829

This is the Infant sampling line which terminates in a 15mm Male-15mm Female straight connector and has a reduced 5mm central bore. This is suitable for patients of between 1kg and 5kg.

Price: £12.50 each



Order code 3827B

This is similar to the original Nomoline and terminates in a single male luer connection. This sampling line is suitable for use in patients weighing less than 1kg and greater than 100kg when used with our connectors and Y-pieces.

Price: £14.50 each



Please note that the Nomoline adaptor (108220) with 15cm line and a single female luer connector is still available.

Price £14.50 each.

All items are for veterinary use only