



WHY BUY A CARDIOSTORE?

BENEFITS

- **Simplicity of use**
- **Mobile and fits in the palm of your hand**
- **Allows user to produce fast professional ECG reports**
- **Saves Time**
 - Allows the ECG's to be faxed or e-mailed in a Windows format to a cardiologist
- **Saves patient's ECG to a ram card uploadable to a Computer for safe keeping**
- **Also monitors heart rate, respiration rate and apnoea.**
- **One ECG machine for a range of patient sizes**




SEE WHAT OUR CUSTOMERS THINK!

"The CardioStore is a very useful portable ECG monitor which is easy to use in the field". **Sue Roberts BVMS CertVC MRCVS - SLR Cardiology Referrals, West Yorkshire.**

SUMMARY

The CardioStore is a six limb lead, 12-bit digital diagnostic recording system with Windows software developed especially for the veterinary market. The CardioStore is hand-held and can move where the patient does, therefore it is useful in theatres due to its size and ability to also measure and notify for signs of apnoea and respiratory rates. The CardioStore allows you to record, pause and rewind whilst the ECG is still taking place. Its Windows software will allow you to store, measure and edit the ECG's for quick diagnosis and reports that can be sent away to a cardiologist at the click of a button. Saving time and space the CardioStore is an amazing piece of monitoring equipment with an incredible price.



Contact us for further information: Tel: 01626 365505 Fax: 0870 129 4705

E-mail: enquiries@vetronic.co.uk

www.vetronic.co.uk

Quality Anaesthetic monitoring
and Ventilation equipment



If you are interested in this product and seek further information,
let us do the work. Write down times you are free and we will do
our best to get back to you as soon as possible at a convenient
time for you.

CARDIOSTORE FAX BACK FORM

Practice Name:	
Person/ Persons interested:	
Address:	
Telephone Number:	
Fax:	
Email:	
Best time for us to contact you regarding your enquiry:	

Vetronic Services Ltd

TEL: 01626 365505

FAX: 0870129 4705

E-mail: enquiries@vetronic.co.uk

www.vetronic.co.uk

Registered Office: 12 Henley's Business Park, Manor Road, Abbotskerswell, Newton Abbot, Devon TQ12 5NF
CRN:4694235; VAT No 585 6341 12. Managing Director: Keith Simpson BVSc MRCVS MIET(Electronics)

*Designed
by Vets, for Vets*