



WHY BUY THE SAPPHIRE PLUS ECG?

BENEFITS

- **Can be used on smaller animals**
 - e.g. Cats
 - **Can be used in dry or water treadmill investigations**
 - **Suitable for long term monitoring. ***
 - **Saves space**
 - Weighing in at just 50 grams and dimensions of 75 x 45 x 18 mm.
 - Saves filling away ECG paper records.
 - **Immediate display of ECG's on a remote computer screen.**
 - **Integrates with the Vetronic Services VitalMonitor Software Suite.**
 - Full wireless capabilities link to PC.
 - Can be added as a module to VitalStore Multimonitor.
 - **Cost Effective**
 - Can print the required sections of the ECG rather than the whole thing.
 - **Clear digital representation**
 - Reduces smart paper reports for presentation to clients.
- *(dependant on battery type)



SEE WHAT OUR CUSTOMERS THINK!

"The Sapphire Plus is a small, compact reliable and versatile ECG system that is affordable for practices and really reliable. I find it a very useful piece of kit as I like the idea of having a telemetry system within the practice and the Sapphire Plus delivers this". **David Sewell BVM&S, BSc MRCVS – Heath Side Veterinary Surgery, Hampshire.**

"This product is really good because unlike other bulky ECG machines with yards of printer paper to analyse, the Sapphire Plus wirelessly transmits the patients data via Bluetooth to the computer system, storing it there for later use. This product is also good for patients who need longer term care or overnight ECG monitoring as it can be attached to the patient and will constantly take readings over longer time frames, storing the data directly on to any Windows compatible computer via Bluetooth." **Rachel Dyer RVN Dip Advanced Vet Nursing (Medicine) – Anderson Moores Veterinary Specialists, Hampshire.**





SUMMARY

Sapphire Plus ECG system is the smallest 6-lead wireless ECG system on the veterinary market. It is a six lead ECG machine weighing in at a mere 50 grams. The patient would almost never know it was there. This product is designed for customers looking to upgrade from a large bulky ECG machine to save space and time. This product can be added as a module to the VitalStore MultiMonitor Software Suite. At 50 grams and with wireless capabilities we are offering a fantastic product at an amazing price.

If you are interested in this product and would like further information, **let us do the work.** Write down the times when you are free and we will do our best to get back to you as soon as possible and at a convenient time for you.

SAPPHIRE PLUS ECG FAX BACK FORM

Practice Name:	
Person/ Persons interested:	
Address:	
Telephone Number:	
Fax:	
Email:	
Best time for us to contact you regarding your enquiry:	

Vetronic Services Ltd. TEL: 01626 365505

FAX: 0870 129 4705 E-mail: enquiries@vetronic.co.uk