



DRY PTFE REQUIREMENTS FOR TAFONIUS SERVICE



For a Tafonius service of any model a full can of DRY PTFE lubricant is required. This Dry lubricant is used to coat the piston diaphragms and improve movement and prevent sticking.

With airline regulations getting more and more stringent it is no longer possible to bring or ship this item ahead and so we need you, the customer, to acquire some of this prior to our visit. Please note that ***WITHOUT THIS PRODUCT WE CANNOT PERFORM THE TAFONIUS SERVICE***

There are a number of manufacturers that produce a suitable Dry PTFE spray. The ones that are currently known to us and have proved suitable are:

CRC Dry lube, in a green can - **our preferred product**

Rocol Dry PTFE Part code: 34235.

Wurth Dry PTFE: Part number is Wurth 0893550

The picture below shows both the Rocol and the CRC products.



The CRC product seems to be popular in Sweden and Norway

Here is a picture of the Wurth product that I have also used in Sweden. This is not as good as the CRC or Rocol products.

